Montavit

Cathejell Lidocain C

Sterile lubricant gel

Product information leaflet

en

Sterilni lubrikant u vidu gela

Letak s informacijama o proizvodu

bs

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





Cathejell Lidocain C

sterile lubricant gel

Cathejell Lidocain C is a sterile, water-soluble, visually clear gel, which is used as a lubricant for the insertion of catheters, endoscopes, probes, or other medical instruments in the urogenital tract to protect the mucosa.

Cathejell Lidocain C is applied topically and acts topically for lubrication purposes, thus facilitating insertion of medical instruments during procedures in the urogenital tract via intraurethral and intravaginal pathways.

Cathejell Lidocain C serves as a lubricant for use on mucous membranes. In addition to its principal intended action, lidocaine hydrochloride is added as a local anaesthetic agent to relieve pain during diagnostic and therapeutic procedures. The anaesthetic effect starts 5 – 10 minutes after application and lasts for 20 – 30 minutes.

To prevent infections within the area of application, Cathejell Lidocain C contains chlorhexidine dihydrochloride. This disinfecting agent acts against a large number of bacteria and some fungi and viruses.

Cathejell Lidocain C is instilled into the urethra prior to the insertion of catheters and other instruments. The gel passes through the urethra ahead of the instrument, gently unfolding and dilating the urethra which would otherwise be collapsed and simultaneously ensuring the formation of a lubricant film between the urethral mucosa and the to-be-inserted instrument.

Cathejell Lidocain C thereby reduces the likelihood of local injuries and the subsequent development of strictures. Decreased friction and optimal gliding performance are accompanied by the reduction of pain and discomfort for the patient.

Due to the above-described preparation of the urethra the patient relaxes the pelvic floor muscles, thereby allowing the instrument to pass more easily through the urethra into the bladder.

Directions for use:

Use in the urethra before inserting a catheter, endoscope, or other medical instruments

Important: Slowly apply (instill) into the urethra before inserting any instruments.

The accordion syringes contain 12.5 g or 8.5 g gel, of which approximately 10 g and 6 g, respectively, are inserted into the urethra during instillation.

Method of administration:

- Clean or disinfect, if appropriate, the external opening of the urethra
- Peel off the paper from the transparent cover up to the tapered end of the blister (Picture 1)
- Snap off the tip with a short sharp jabbing action into the blister pack (Picture 2)
- Remove the tip completely so that it cannot be inserted accidentally into the urethra
- Squeeze out one drop of gel so that the application nozzle can be inserted more easily (Picture 3)
- Slowly insert the gel by applying moderate pressure on the syringe. (Picture 4)
- Empty the accordion syringe completely and remove it pressed (Picture 5)

- The syringes are intended for single use; discard any remaining gel.

Use with endoscopes and probes

The accordion syringes contain 12.5 g or 8.5 g gel, of which the designated amount of up to about 10 g and 6 g, respectively, is applied.

- Peel off the paper from the transparent cover up to the tapered end of the blister (Picture 1)
- Snap off the tip with a short sharp jabbing action into the blister pack (Picture 2)
- Remove the tip completely so that it cannot be inserted by mistake
- Slowly apply the gel by exerting moderate pressure on the syringe.

Apply gel to the probe / endoscope / instrument depending on the intended use and spread evenly.

To prevent drying out, the gel should not be applied to the instrument until immediately before use.

The syringes are intended for single use; discard any remaining gel.

Dosage instructions:

The syringes contain 12.5 g or 8.5 g gel, of which about 10 g and 6 g, respectively, can be removed.

The amount of gel to be used depends on the patient's anatomical characteristics and the type of endoscope, catheter, tube, probe or instrument inserted.

For catheterisation, never instill more than 1 syringe at any one time. The contents of one syringe are sufficient to fill the urethra.

The maximum dose is 16g Cathejell Lidocain C – anaesthetising gel for adults of normal body weight.

Cathejell Lidocain C is indicated in adults, adolescents and children aged 2 to 11 years.

In children below 2 years of age, Cathejell Lidocain C must not be used.

In children below 12 years of age, doses must not exceed 2.9 mg lidocaine (≈ 0.15 g gel or 0.14 ml gel) / kg BW. No more than four doses shall be administered within 24 hours.

In patients at risk, elderly, debilitated and acutely ill patients and in cases of impaired liver function or severe kidney dysfunction, inflammation or ulcer at the administration site, the dosage should be adjusted accordingly. The maximum dose must be calculated in mg lidocaine hydrochloride/kg body weight (2.9 mg lidocaine hydrochloride/kg BW).

Do not exceed the recommended dose.

Do not use Cathejell Lidocain C

- If the patient is allergic (hypersensitive) to the active substances or any of the other ingredients of Cathejell Lidocain C,
- If the patient is allergic (hypersensitive) to certain other local anaesthetics (of the amide type),
- in children below 2 years of age,
- in bulbocavernosus reflux (an injury to the thin urethral lining, which may cause the lubricant to infiltrate into erectile tissue, where absorption may occur),
- in patients with a severely weak heart, extremely slow heartbeat, heart conduction disturbances (AV block), shock due to heart failure or reduced blood volume.

Take special care with Cathejell Lidocain C

- If the patient has an inflammation, injury, or an ulcer of the local epithelium where Cathejell Lidocain C is administered
- If the patient's liver or kidney function is severely impaired
- If the patient has heart or respiratory tract dysfunction
- If the patient is elderly, in a weak state of health or acutely ill
- If the patient is prone to seizures or epilepsy
- If the patient is suffering from a certain muscle disease (myasthenia gravis)
- If the patient is on treatment with certain medicines for heart rhythm disorders, known as class III antiarrhythmics (e.g. amiodarone), as the effects on the heart may be enhanced
- in children, as absorption of lidocaine into the blood stream may be increased.
- in patients with acute porphyria or methemoglobinemia. In this case, medical advice should be sought before use.

If more than one syringe is instilled or if the local epithelium is ulcerous / inflamed, absorption of lidocaine may be increased. This may lead to an overdose, particularly in children and elderly patients, with subsequent disorders of the central nervous and cardiovascular system (see also "side effects").

Cathejell Lidocain C must be applied by a doctor or an experienced medical staff.

Cathejell Lidocain C must not come into contact with the eyes, middle ear or wounds.

In general anaesthesia, a lubricant without lidocaine should be preferred.

Interaction between Cathejell Lidocain C and other substances:

The use of anionic tensides (certain soaps) can considerably reduce the disinfectant action of chlorhexidine.

Due to possible increased effects on the heart, Cathejell Lidocain C should not be used at the same time as:

Medicines containing lidocaine or other amide-type local anaesthetics, antiarrhythmics (used to treat heart rhythm disorders) or beta-blockers or calcium channel blockers (used for high blood pressure).

Taking cimetidine (inhibits stomach acid production) can increase the risk of side effects.

Pregnancy and breast-feeding:

During pregnancy and breast-feeding, Cathejell Lidocain C should be used only after careful consideration of the benefits and risks by a doctor and after the individual dose has been established.

Repeated use during pregnancy and breast-feeding is not recommended.

Please note that there should be a 12-hour interval between using Cathejell Lidocain C and subsequent breast-feeding.

Women planning a pregnancy or who think they are pregnant should tell their doctor.

Driving and using machines

Effects on the ability to drive and use machines are unlikely but cannot be completely ruled out in cases of increased individual sensitivity.

Side effects:

Side effects rarely occur after application of Cathejell Lidocain C, provided that the product is used according

to the dosage recommendations/recommendations for use and the necessary precautions are taken.

In rare cases, allergic reactions such as redness, burning or itching may occur. However, more severe, general hypersensitivity reactions and even shock are also possible. In view of the low exposure levels and short duration of use with Cathejell Lidocain C, other significant side effects/interactions are unlikely.

Overdose:

An overdose as well as high plasma levels due to rapid absorption of lidocaine hydrochloride or reduced tolerance can cause systemic side effects. The following symptoms may occur:

- CNS: nervousness, dizziness, blurred vision, tremor. These signs may not necessarily occur. In some patients, intoxication manifests as drowsiness, unconsciousness, or respiratory arrest.
- Cardiovascular: hypotension, bradycardia, asystole.

If signs of an overdose appear/emerge, consult a doctor immediately. The treatment of an overdose is symptomatic, e.g. by administering anticonvulsants, and/or performing supporting cardiopulmonary emergency measures.

Cathejell Lidocain C may cause methemoglobinemia.

Special warnings for safe use:

- In anaesthesia, a lubricant without lidocaine should be preferred.
- Only use when sterile paper is intact and when gel is clear.
- Only use once. If re-used or re-sterilized, functionality and product quality can not be guaranteed.
- After instillation, a 5 to 10-minute interval should be respected before pain-free insertion of the instrument is possible.

Composition:

1 g contains 20 mg lidocaine hydrochloride and 0.5 mg chlorhexidine dihydrochloride.

The other ingredients are: hydroxyethylcellulose, glycerol, sodium lactate, water for injections.

Disposal:

After the first use, the syringe must be considered contaminated and should be disposed of in accordance with national regulations.

Pack sizes:

1 x 12.5 g, 5 x 12.5 g, 25 x 12.5 g

1 x 8.5 g, 5 x 8.5 g, 25 x 8.5 g

Shelf life:

As printed on the pack.

Storage instructions:

Do not store above 25°C.

Keep the blister in the carton to protect from light.

Keep out of reach of children.

Keep dry.

Reporting of incidents

In case of serious incidents during the use of Cathejell please inform the national competent authority and the manufacturer.

A summary of safety and clinical performance can be found on the European database on medical devices:

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